


Silicone rubber implants - implants for episcleral surgery - sterile disposable

PRODUCT LINE AND CODES			
CODE	REF	SIZE	SHAPE
RUBBER STRIPES (STRIPE)			
BPP1 *	MMD-701	2x0,8 mm	0.8 
BP1 *	MMD-703	2,5x0,7 mm	0.7 
BM1 *	MMD-705	3,5x0,9 mm	0.9 
BL1 *	MMD-707	4,2x1,1 mm	1.1 

* Forniti di sleeve

PRODUCT LINE AND CODES			
CODE	REF	SIZE	SHAPE
RUBBER STRIPES (TYRE)			
TBM6	MMD-736	2,5x7 mm	
TBM7	MMD-737	2,5x7 mm	
TBM9	MMD-739	2,5x9 mm	
MACULAR BUCKLE			
BMM4	MMD-726	5 mm	

TECHNICAL DATA

	Rubber Stripes (Stripe)	Macular Buckle and Rubber Stripes (Tyre)
<i>Silicone Rubber</i>	XIAMETER RBB-2100-30/70	PTMED45 Clear

KIND OF PACKAGING

Envelope content sterilized at Ethylene Oxide
Small Carton Box (90x130x35)mm